

# A short introduction in the EUROIMMUN program (www.ifq-portal.de) 1<sup>st</sup> step: User Registration

IfQ-Lübeck

Institut für  
Qualitätssicherung  
Lübeck



An institution of the EUROIMMUN Medizinische Labordiagnostika AG


## Quality Assessment Portal

### Information

- [Contact / site notice](#)
- [Flyer](#)
- [Conditions of participation](#)
- [Time table and price list](#)
- [Reference laboratories](#)
- [Product catalogue](#)
- [Feedback form](#)
- [English software manual](#)
- [软件使用说明 \(中文\)](#)
- [Accreditation document](#)




### Login

Language: English  ▼

E-Mail

Password

Login

Not registered? 

Password forgotten?

New registration



# 1<sup>st</sup> step: User Registration

## User registration

### User profile

Please insert your data and click the "Register"-button, if you want do participate on the EUROIMMUN AG quality assessments. Your E-Mail address will be your username(login name) for the usage of the Quality-Assessment-Portal.

Please don't register multiple times. With this registration you accept the conditions of participation!

<b>Firstname:</b>	<input type="text"/>	<b>Salutation:</b>	<input type="text" value="generally"/>
<b>Surname:</b>	<input type="text"/>	<b>Title:</b>	<input type="text" value="no title"/>
<b>Street:</b>	<input type="text"/>	<b>Institution:</b>	<input type="text"/>
<b>Postcode &amp; city:</b>	<input type="text"/>	<b>Department:</b>	<input type="text"/>
<b>Country:</b>	<input type="text" value="Choose country"/>	<b>Password:</b>	<input type="text"/>
<b>Tel. number:</b>	<input type="text"/>	<b>Password(rep.):</b>	<input type="text"/>
<b>Fax number:</b>	<input type="text"/>		
<small>(optional)</small>			
<b>E-Mail:</b>	<input type="text"/>		
<b>E-Mail (rep.):</b>	<input type="text"/>		

### Billing address

If bill address is the same as the above address you can easily copy it from above.

Copy address

<b>Firstname:</b>	<input type="text"/>	<b>Postcode &amp; city:</b>	<input type="text"/>
<b>Surname:</b>	<input type="text"/>	<b>Department:</b>	<input type="text"/>
<b>Street:</b>	<input type="text"/>	<b>Tax number:</b>	<input type="text"/>
		<small>(optional)</small>	

Register

Cancel

Please enter address, to which sera should be sent to and name which should be listed on the certificate.

Click on Copy address, if shipment address and invoice address is the same

Click register to close. Password will be sent to registered Email address

# Registration for a special scheme

Please login with your Email address and password

If you would like to change your password or address you can click here

Logout here

Read the description of the scheme

The following quality assessments are available

Name	Registration status				
	I/2017		II/2017		
Borrelia <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
CCP <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
ANCA <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Skin <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Arboviruses <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Thyroid gland <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Coeliac disease <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Phospholipids <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
ANA I <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
EBV <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Allergology <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
HSV <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Parvovirus <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
ANA II <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Liver <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
kidney <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
HEV <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>
Neuronal <a href="#">Information</a>	Not registered	[O]	Not registered	[O]	<a href="#">Register/Deregister</a>

[O] -> open for registration [C] -> registration closed [S] -> shipped [E] -> in evaluation



# Registration for a quality assessment

Thyroid gland	Information	Not registered [0]	Not registered [0]	Register/Deregister
<b>Offered sets of specimen:</b>	1 set of specimen (3x300µl) 2 sets specimen (3x600µl)			
<b>Content and purpose:</b>	This quality assessment program shall contribute to quality control and standardization of the determination of autoantibodies against thyroid gland of your diagnostic laboratory. The quality assessment includes the evaluation of antibodies against TPO, Tg and the determination of TRAb. The test systems are evaluated individually. A certificate is only awarded upon correct analysis of all three samples.			

# Registration for a quality assessment

## Register/Deregister

This quality assessment will be charged. You can get the specific costs for the different sets of specimen from your local distributor.

Time schedule

Find distributor

Information about the quantity of the different sets of specimen for the quality assessment can be found in the time schedule.

I accept the conditions of participation for all quality assessment participations of the Institute for Quality Assurance Lübeck. ([Conditions of participation](#))

### Thyroid gland I/2017:

Register until: Nov 30, 2016

Shipping date: May 2, 2017

End date: May 30, 2017

Report date: Jun 27, 2017

Choose set of specimen!  1 set of specimen  2 sets of specimen

Order number: (optional) \*

Register (with costs)

### Thyroid gland II/2017:

Register until: Jun 30, 2017

Shipping date: Oct 17, 2017

End date: Nov 14, 2017

Report date: Dec 12, 2017

Choose set of specimen!  1 set of specimen  2 sets of specimen

Order number: (optional) \*

Register (with costs)

\*If you want to use an order number for this quality assessment, you can enter it here.

Close

Click on 1 set or 2 sets of specimen depending on how many volume you need for your testings

You only have to fill in the order number if you need it

Click on register

# Registration for a quality assessment

## Register/Deregister

This quality assessment will be charged. You can get the specific costs for the different sets of specimen from your local distributor.

Time schedule

Find distributor

Information about the quantity of the different sets of specimen for the quality assessment can be found in the time schedule.

### Thyroid gland I/2017:

Shipping date: May 2, 2017  
End date: May 30, 2017  
Report date: Jun 27, 2017  
Status: **registered** 1 set of specimen (3x300µl)  
Order number:

Deregister

### Thyroid gland II/2017:

Register until: Jun 30, 2017  
Shipping date: Oct 17, 2017  
End date: Nov 14, 2017  
Report date: Dec 12, 2017  
Choose set of specimen!  1 set of specimen  2 sets of specimen  
Order number: (optional) \*

Register (with costs)

\*If you want to use an order number for this quality assessment, you can enter it here.

Close

# Registration for a quality assessment

## The following quality assessments are available

Name		Registration status		
		I/2017	II/2017	
Borrelia	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
CCP	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
ANCA	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Skin	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Arboviruses	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Thyroid gland	<a href="#">Information</a>	registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Coeliac disease	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Phospholipids	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
ANA I	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
EBV	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Allergology	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
HSV	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Parvovirus	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
ANA II	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Liver	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
kidney	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
HEV	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>
Neuronal	<a href="#">Information</a>	Not registered [O]	Not registered [O]	<a href="#">Register/Deregister</a>

[O] -> open for registration [C] -> registration closed [S] -> shipped [E] -> in evaluation

Your registration for the first half year for thyroid gland was successful

# Data entry (after shipment of the sera)

Home   Quality assessments   User profile   Site notice   Contact

Quality assessment registration   Enter methods/results   Evaluations   Refresh   Change language

Home

### New informations concerning our quality assessment

September 20, 2016:  
[Letter\\_QA\\_2017.pdf](#)

March 1, 2016:  
Ladies and gentlemen,  
Please note the following changes:

### The following quality assessments are available

Name		Registration status		
		I/2017	II/2017	
Thyroid gland	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>
Neuronal	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>
ANCA	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>
kidney	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>
Liver	<input type="button" value="Information"/>	[1x] <b>registered</b> [S]	Not registered [O]	<input type="button" value="Register/Deregister"/>
Skin	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>
CCP	<input type="button" value="Information"/>	Not registered [C]	Not registered [O]	<input type="button" value="Register/Deregister"/>

„Enter methods/results“  
for data entry





# Data entry

## Quality assessments

Quality assessment registration   Running quality assessment   Completed

Name	Shipping date	Registration date	Number of Assays	Number of results	
Thyroid gland, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry
Phospholipids, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry
ANCA, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry
Coeliac disease, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry
CCP, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry
Skin, May 2014	May 13, 2014	Oct 25, 2013	0	0 / 0	Information Data entry

Number of defined assays

Number of entered results

Click on „Data entry“

# 3<sup>rd</sup> step: Data entry

## Define methods

Quality assessments » Edit/add methods for "Phospholipids, May 2014"

### Quality assessment Information

<b>Name:</b>	Phospholipids, May 2014	<b>Phase:</b>	result input
<b>Register until:</b>	Nov 30, 2013	<b>Participant status</b>	registered
<b>Shipping date:</b>	May 13, 2014	<b>Registration date:</b>	Oct 25, 2013
<b>End date:</b>	Jun 10, 2014	<b>Set of specimen:</b>	2 sets (3x300µl)
		<b>Order number:</b>	

**Content and purpose:** This quality assessment program shall contribute to quality control and standardization of the determination of autoantibodies against phospholipids of your diagnostic laboratory. The quality assessment includes the evaluation of antibodies against Cardiolipin and  $\beta$ 2-Glycoprotein.

### Quality assessment Edit/add methods for "Phospholipids, May 2014"



Add assay from past participation.



Add assay

You can click here to copy  
your assays from the last QA

Please click here  
to define new  
assays

## 3<sup>rd</sup> step: Data entry – Define methods

**Add assay**

**Select already used assay:**  
Choose assay

or add new assay:

**Format:** Choose

**Manufacturer:** Choose

+ Add assay    Cancel

Enter result    - Remove all assays

Please choose method  
and manufacturer

## 3<sup>rd</sup> step: Data entry- Define methods

**Add assay**

Select already used assay:

or add new assay:

**Format:**

**Manufacturer:**

**Parameter:**  
 anti-cardiolipin  
 anti-β2-glycoprotein

**Conjugates:**  IgG  IgM  IgAGM

**Name of manufacturer:**

**Name of assay:**

Add assay

## 3<sup>rd</sup> step: Data entry

Quality assessment Edit/add methods for "Phospholipids, May 2014"

<b>Format:</b>	ELISA	<input type="button" value="- Remove"/>
<b>Manufacturer:</b>	xyz	<input type="button" value="📄 Replace"/>
<b>Name:</b>	Anti-Cardiolipin-ELISA	
<b>Parameter:</b>	anti-cardiolipin,	
<b>Conjugates:</b>	IgG	

Click here to enter results

Add next assay

# Adding or changing of methods

Save Edit/add methods Refresh Download result preview Back

Quality assessments » Results for "Phospholipids, May 2014"

## Results for "Phospholipids, May 2014"

Sample	Name of assay	Parameter	Qual. result	Cut-Off ?	Quant. result	Unit	perform in-house ?
Sample 1	Anti-Cardiolipin xyz (IgG)	anti-cardiolipin	Choose			Choose	<input type="radio"/> Yes <input type="radio"/> No
Sample 2	Anti-Cardiolipin xyz (IgG)	anti-cardiolipin	Choose			Choose	<input type="radio"/> Yes <input type="radio"/> No
Sample 3	Anti-Cardiolipin xyz (IgG)	anti-cardiolipin	Choose			Choose	<input type="radio"/> Yes <input type="radio"/> No

You can add or change existing methods whenever you want to.

# Save data

Save Edit/add methods Refresh Download result preview Back

Quality assessments » Results for "Phospholipids, May 2014"

**Results for "Phospholipids, May 2014"**

Sample	Name of assay	Parameter	Qual. result	Cut-Off ?	Quant. result	Unit	perform <sup>?</sup> in-house
Sample 1	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	positive	12	120	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 2	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	negative	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 3	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	negative	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No

Click on „Save“ after data entry to save your results.

(recommended after each sample, because it also causes the take over of Cut-off and units)

# Result preview

Save Edit/add methods Refresh Download result preview Back

Quality assessments » Results for "Phospholipids, May 2014"

### Results for "Phospholipids, May 2014"

Sample	Name of assay	Parameter	Qual. result	Cut-Off <sup>?</sup>	Quant. result	Unit	perform <sup>?</sup> in-house
Sample 1	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	positive	12	>120	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 2	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	negative	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 3	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	negative	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No

Results are saved

Click here for an overview of your results

You can also enter > or < if you are not diluting



# Download result preview

## Quality assessment QV 1632-0303-1, May 2014, "Phospholipids": Preview of your results

Sample	Name of assay	Parameter	Cut-Off	Result		Unit	performed in-house
				qualitative	quantitative		
Sample 1	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	12	positive	>120	GPLU	Yes
Sample 2	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	12	negative	2	GPLU	Yes
Sample 3	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	12	negative	2	GPLU	Yes

You can check your results and save them

# Data entry of other running schemes

Save

Edit/add methods

Refresh

Download result  
preview

Back


Quality assessments » Results for "Phospholipids, May 2014"

## Results for "Phospholipids, May 2014"

Sample	Name of assay	Parameter	Qual. result	Cut-Off <sup>?</sup>	Quant. result	Unit	perform <sup>?</sup> in-house
Sample 1	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	<input type="text" value="positive"/>	12	>120	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 2	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	<input type="text" value="negative"/>	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No
Sample 3	Anti-Cardiolipin-ELISA, ELISA xyz (IgG)	anti-cardiolipin	<input type="text" value="negative"/>	12	2	GPLU	<input checked="" type="radio"/> Yes <input type="radio"/> No

Go back to Quality assessment, if you want to enter results for other schemes

# Data entry of other running schemes

Home Quality assessments User profile Site note Contact 

Refresh

### Quality assessments

Quality assessment registration Running quality assessment Evaluations

Name	Shipping date	Registration date	Number of Assays	Number of results		
EBV, April 2014	Mar 25, 2014	Oct 25, 2013	2	2 / 18	Information	Data entry
Liver, April 2014	Mar 25, 2014	Oct 25, 2013	6	9 / 27	Information	Data entry
ANA, March 2013	Mar 19, 2013	Oct 2, 2012	3	0 / 15	Information	Data entry
Borrelia, March 2013	Mar 19, 2013	Oct 2, 2012	4	12 / 12	Information	Data entry
HSV, March 2013	Mar 19, 2013	Oct 2, 2012	0	0 / 0	Information	Data entry
Parvovirus, March 2013	Mar 19, 2013	Oct 2, 2012	0	0 / 0	Information	Data entry

6 entries       (1 of 1)

For this scheme all results for the defined methods are entered

No methods defined

Methods defined, no results are entered



# Access evaluations

When evaluation is ready you will be informed by Email. You can then login again. Click on evaluations

Quality assessments

Quality assessment registration   Running quality assessment   **Evaluations**

Name	Correct assay number	Status	
ANCA, October 2013	5 / 5	certificate	Show evaluation
Phospholipids, October 2013	4 / 4	certificate	Show evaluation
ANA, September 2013	20 / 20	certificate	Show evaluation

3 entries   10   20   50   100   (1 of 1)

After that click on show evaluation

# Evaluations

Quality assessments » Evaluation for "ANA, September 2013"

Evaluation - "ANA, September 2013" [Show information](#)

[Evaluation report](#) [Images](#) [Reference Value](#)

Test systems

▶ Antibodies against cell nuclei (IgG), EUROIMMUN, IFA	Certificate: Yes 3 / 3 correct results
▶ Anti-dsDNA-NcX ELISA , EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-Nucleosomes ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-Centromeres ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-SS-A ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-SS-B ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-Sm ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-nRNP/Sm ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-Scl-70 ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-Jo1 ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Anti-PM-Scl ELISA (IgG), EUROIMMUN, ELISA	Certificate: Yes 3 / 3 correct results
▶ Crithidia lucilia (dsDNA), EUROIMMUN, IFA	Certificate: Yes 3 / 3 correct results

Evaluation report contains clinical details of samples, rates of correct results, numbers of participants and an overview of test systems used

IFA pictures or EuroLineScan pictures of QA samples

Results of reference laboratories

Still questions?

If there are still questions mail to  
[ifq@euroimmun.de](mailto:ifq@euroimmun.de)