



Time Schedule 2019

Autoimmunity

Scheme	Format	Dates QAS 2019/I	Dates QAS 2019/II
Autoantibodies against cell nuclei Cell nuclei, centromeres, CenpA, CenpB, dsDNA, nucleosomes, ribosomal P protein, RNP/Sm, RNP, Scl-70, Sm, SS-A, SS-B	1 set of specimen 2x400µl	<u>Registration until:</u> November 30, 2018	<u>Registration until:</u> June 30, 2019
Antibodies in autoimmune Myopathies cN-1A, EJ, Ku, Jo-1, Mi-2, Mi-2α, Mi-2β, OJ, PM-Scl, PL-7, PL-12, SRP	1 set of specimen 2 sets of specimen 2x200µl 2x400µl	<u>Shipment:</u> March 19, 2019	<u>Shipment:</u> September 17, 2019
Antibodies in autoimmune liver diseases AMA, ASMA, F-actin, gp210, LC-1, LKM-1, M2, nuclear dots, nuclear membrane, SLA/LP, Sp100	1 set of specimen 2 sets of specimen 2x200µl 2x400µl	<u>Deadline:</u> April 16, 2019	<u>Deadline:</u> October 15, 2019
Autoantibodies against thyroid gland TG, TPO, TRAb	1 set of specimen 2 sets of specimen 2x300µl 2x600µl		
Autoantibodies against neuronal antigens Amphiphysin, Aquaporin 4, CASPR2, CV2, GAD, Hu, LGI1, Ma/Ta, MOG, NMDA-Rezeptor, Recoverin, Ri, SOX1, Titin, Yo, Zic4	1 set of specimen 2 sets of specimen 2x200µl 2x400µl		
Autoantibodies against granulocytes ANCA, MPO, PR3	1 set of specimen 2 sets of specimen 2x200µl 2x400µl		
Autoantibodies against structural proteins of the skin BP 180, BP 230, desmoglein 1 und 3, collagen type VII, desmosomes, epidermal basement membrane	1 set of specimen 2 sets of specimen 2x200µl 2x400µl	<u>Registration until:</u> November 30, 2018	<u>Registration until:</u> June 30, 2019
Antibodies against CCP	1 set of specimen 2 sets of specimen 2x200µl 2x400µl	<u>Shipment:</u> April 30, 2019	<u>Shipment:</u> October 15, 2019
Autoantibodies against phospholipids Cardiolipin (IgG, IgM, IgAGM); β2-Glykoprotein (IgG, IgM, IgAGM)	1 set of specimen 2 sets of specimen 2x150µl 2x300µl	<u>Deadline:</u> May 28, 2019	<u>Deadline:</u> November 12, 2019
Antibodies against tissue transglutaminase (endomysium) and gliadin (IgA, IgG)	1 set of specimen 2 sets of specimen 2x200µl 2x400µl		
Antibodies in autoimmune kidney diseases GBM, PLA2R, THSD7A	1 set of specimen 2 sets of specimen 2x200µl 2x400µl		
Autoantibodies in pernicious anaemia Intrinsic factor, parietal cells	1 set of specimen 2 sets of specimen 2x200µl 2x400µl		



Infectious serology

Scheme	Format	Dates QAS 2019/I	Dates QAS 2019/II
Antibodies against Borrelia burgdorferi sensu lato (IgG, IgM, serological diagnosis)	1 set of specimen 2x250µl 2 sets of specimen 2x500µl	<u>Registration until:</u> November 30, 2018 <u>Shipment:</u> March 19, 2019 <u>Deadline:</u> April 16, 2019	<u>Registration until:</u> June 30, 2019 <u>Shipment:</u> September 17, 2019 <u>Deadline:</u> October 15, 2019
Antibodies against Herpes-simplex-Virus HSV 1, HSV 2, HSV 1/2 (IgG, IgM)	1 set of specimen 2x200µl 2 sets of specimen 2x400µl		
Antibodies against Parvovirus-B19 (IgG, IgM)	1 set of specimen 2x200µl 2 sets of specimen 2x400µl		
Antibodies against Epstein-Barr-Virus EBV-CA (IgG, IgM), EBNA-1 (IgG), Avidity, serological diagnosis	1 set of specimen 2x200µl 2 sets of specimen 2x400µl		
Antibodies against Hepatitis-E-Virus (IgG, IgM, IgAGM)	1 set of specimen 2x200µl 2 sets of specimen 2x400µl		
Antibodies against Zika virus (IgG, IgM, IgA, IgAM)	1 set of specimen 2x300µl 2 sets of specimen 2x600µl		
Antibodies against Dengue virus (IgG, IgM)	1 set of specimen 2x300µl 2 sets of specimen 2x600µl		
Antibodies against Chikungunya virus (IgG, IgM)	1 set of specimen 2x300µl 2 sets of specimen 2x600µl		

Allergology

Scheme	Format	Dates QAS 2019/I	Dates QAS 2019/II
Allergology total IgE	1 set of specimen 2x500µl 2 sets of specimen 2x1000µl	<u>Registration until:</u> November 30, 2018 <u>Shipment:</u> March 19, 2019 <u>Deadline:</u> April 16, 2019	<u>Registration until:</u> June 30, 2019 <u>Shipment:</u> September 17, 2019 <u>Deadline:</u> October 15, 2019
Allergology specific IgE	1 set of specimen 2x1000µl		